



IHE Integration Statement

Rialto 2.1

March 1, 2010

Karos Health Incorporated
295 Hagey Boulevard
1st Floor, West Entrance
Waterloo, ON
N2L 6R5
Canada

IHE Integration Statement

Vendor	Product Name	Version	Date
Karos Health Incorporated	Rialto	2.1	1 March 2010

This product implements all of the transactions specified in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options specified below:

Integration Profiles Implemented	Actors Implemented	Options Implemented
Access to Radiology Information (ARI)	Image Manager/Image Archive Report Repository	none none
Antepartum Education (APE)	Content Consumer	View Option
Antepartum History and Physical (APHP)	Content Consumer	View Option
Antepartum Laboratory (APL)	Content Consumer	View Option
Antepartum Summary (APS)	Content Consumer	View Option
Audit Trail and Node Authentication (ATNA)	Audit Record Repository Secure Application Secure Node	none Encryption Encryption
Basic Patient Privacy Consents (BPPC)	Content Consumer Document Consumer	Basic Patient Privacy Acknowledgement View, Basic Patient Privacy Enforcement Basic Patient Privacy Proof
Consistent Presentation of Images (CPI)	Image Manager/Image Archive	none
Composite Triage and Nursing Note (CTNN)	Content Consumer	View Option
Consistent Time (CT)	Time Client	none
ED Physician Note (EDPN)	Content Consumer	View Option
Evidence Documents (ED)	Image Manager/Image Archive	none
Emergency Department Referral (EDR)	Content Consumer	View Option
Key Image Note (KIN)	Image Manager/Image Archive	none
Immunization Content (IC)	Content Consumer	View Option
Labor and Delivery Summary (LDS)	Content Consumer	View Option
Labor and Delivery History and Physical (LDHP)	Content Consumer	View Option
Maternal Discharge Summary (MDS)	Content Consumer	View Option
Nursing Note (NN)	Content Consumer	View Option
Nuclear Medicine Image (NM)	Image Manager/Image Archive	none
Patient Demographics Query (PDQ)	Patient Demographics Consumer Patient Demographics Supplier	none none
Patient Demographics Query HL7 V3 (PDQV3)	Patient Demographics Consumer Patient Demographics Supplier	none none

Patient Information Reconciliation (PIR)	ADT Patient Registration Image Manager/Image Archive Performed Procedure Step Manager	none none none
Patient Identifier Cross-referencing for MPI (PIX)	Patient Identifier Cross-reference Consumer Patient Identity Source Patient Identity Cross-reference Manager	PIX Update Notification, none, none
Patient Identifier Cross-Reference HL7 V3 (PIXV3)	Patient Identifier Cross-reference Consumer Patient Identity Source Patient Identity Cross-reference Manager	PIXV3 Update Notification, none, none
Patient Plan Care (PPOC)	Content Consumer	View Option
Personnel White Pages (PWP)	Personnel White Pages Consumer Personnel White Pages Directory	none none
Simple Image and Numeric Report (SINR)	Report Repository	none
Sharing Value Sets (SVS)	Value Set Consumer Value Set Repository	HTTP Binding - Web Services Binding none
Scheduled Workflow (SWF)	ADT Patient Registration Image Archive/ Image Manager Performed Procedure Step Manager	none none none
Triage Note (TN)	Content Consumer	View Option
Sharing of Lab Report (XD-LAB)	Content Consumer	View Option
Cross Community Access (XCA)	Initiating Gateway Responding Gateway	none none
Cross-Enterprise Document Reliable Interchange (XDR)	Document Source	Basic Patient Privacy Enforcement, Multiple Document Submission
Cross-Enterprise Clinical Documents Sharing (XDS.b)	Document Consumer Document Registry Document Repository Document Source Patient Identity Source	Basic Patient Privacy Enforcement, Basic Patient Privacy Proof Patient Identity Feed (HL7 V2) none Document Addendum, Document Replace, Document Transformation, Folder Management Patient Identity Feed (HL7 V2)
Cross-enterprise Document Sharing for Imaging (XDS-I.b)	Document Consumer Document Registry Document Repository Imaging Document Consumer Imaging Document Source	none none none none PDF report, Set of DICOM Instances, Text Report
Cross-Enterprise Sharing of Medical Summaries (XDS-MS)	Content Consumer	View Option

Cross-Enterprise Document Sharing of Scanned Documents (XDS-SD)	Content Consumer	View Option
Exchange of Personal Health Record Content (XPHR)	Content Consumer	View Option
Cross Enterprise User Assertion (XUA)	X-Service Provider X-Service User	none none
Internet address for vendor`s IHE information: http://www.karoshealth.com/standards		
Links to Standards Conformance Statements for the Implementation		
HL7	http://www.karoshealth.com/standards	
DICOM	http://www.karoshealth.com/standards	
Links to general information on IHE		
In North America: http://www.ihe.net	In Europe: http://ihe-europe.org	In Japan: http://www.jira-net.or.jp/ihe-j