

IHE Integration Statement

Vendor	Product Name	Version	Date
Karos Health	Rialto for Diagnostic Imaging	2.0	January 19, 2010

This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below:

Profiles Implemented	Actors Implemented	Options Implemented
Antepartum Education	Content Consumer	View Option
Antepartum History and Physical		
Antepartum Laboratory		
Antepartum Summary		
Audit Trail and Node Authentication	Secure Application	ATNA Encryption
Basic Patient Privacy Consents	Content Consumer	Basic Patient Privacy Acknowledgement View
	Document Consumer	Basic Patient Privacy Enforcement, Basic Patient Privacy Proof
C38 Patient Level Quality Document Using IHE Medical Summary (XDS-MS)	Content Consumer	
Composite Triage and Nursing Note		View Option
Consistent Time	Time Client	
Cross-Enterprise Clinical Documents Share (XDS.b)	Document Consumer	Basic Patient Privacy Enforcement
	Document Source	
Cross-Enterprise Document Reliable Interchange		
Cross-Enterprise Document Sharing for Imaging	Document Consumer	
	Imaging Document Consumer	
	Imaging Document Source	
Cross-Enterprise Document Sharing of Scanned Documents	Content Consumer	View Option
Cross-Enterprise Sharing of Medical Summaries		
Cross Enterprise User Authentication	X-Service User	

Internet address for vendor's IHE information: http://www.karoshealth.com/Rialto_Standards.htm

Links to Standards Conformance Statements for the Implementation

HL7	
DICOM	http://karoshealth.com/dicom_cs_rialto_20.pdf

Links to general information on IHE

In North America: www.ihe.net	In Europe: www.ihe-europe.net	In Japan: www.ihe-j.org
--	---	--

IHE Integration Statement

Vendor	Product Name	Version	Date
Karos Health	Rialto for Diagnostic Imaging	2.0	January 19, 2010

This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below:

Profiles Implemented	Actors Implemented	Options Implemented	
ED Physician Note	Content Consumer	View Option	
Emergency Department Encounter Summary			
Emergency Department Referral			
EMS Transfer of Care			
Exchange of Personal Health Record Content			
HITSP Collect and Communicate Audit Trail	Secure Application		
HITSP Emergency Encounter Summary Documents Component	Content Consumer		
HITSP Encounter Document Using IHE Medical Summary (XDS-MS)			
HITSP Entity Identity Assertion Component	X-Service Provider		
HITSP Immunization Content	Content Consumer		
			Document Consumer
			Document Source
HITSP Patient Demographics Query Transaction	Patient Demographics Consumer		
HITSP Reliable Document Exchange	Document Recipient		
HITSP Secured Communication Channel Transaction	Secure Application		
HITSP Sharing Radiology Results Transaction Package	Document Consumer		
	Document Source		

Internet address for vendor's IHE information: http://www.karoshealth.com/Rialto_Standards.htm

Links to Standards Conformance Statements for the Implementation

HL7	
DICOM	http://karoshealth.com/dicom_cs_rialto_20.pdf

Links to general information on IHE

In North America: www.ihe.net	In Europe: www.ihe-europe.net	In Japan: www.ihe-j.org
--	---	--

IHE Integration Statement

Vendor	Product Name	Version	Date
Karos Health	Rialto for Diagnostic Imaging	2.0	January 19, 2010

This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below:

Profiles Implemented	Actors Implemented	Options Implemented	
HITSP Sharing Radiology Results Transaction Package	Imaging Document Consumer		
	Imaging Document Source		
HITSP T16: Consistent Time	Time Client		
HITSP XDS-SD	Content Consumer		
Immunization Content			View Option
Labor and Delivery History and Physical			
Labor and Delivery Summary			
Lab Report Document Component			
Manage Sharing of Documents Transaction Package with Doc Integrity Check	Document Consumer		
	Document Source		
Manage Sharing of Documents (XDS.B Option)	Document Consumer		
	Document Source		
Maternal Discharge Summary	Content Consumer	View Option	
Nursing Note			
Patient Demographic Query HL7 V3	Patient Demographics Consumer		
Patient Demographics Query			
Patient ID Cross-Referencing using IHE PIX	Patient Identifier Cross-reference Consumer		
Patient Identifier Cross-Reference HL7 V3			
Patient Identifier Cross-referencing for MPI			

Internet address for vendor's IHE information: http://www.karoshealth.com/Rialto_Standards.htm

Links to Standards Conformance Statements for the Implementation

HL7	
DICOM	http://karoshealth.com/dicom_cs_rialto_20.pdf

Links to general information on IHE

In North America: www.ihe.net	In Europe: www.ihe-europe.net	In Japan: www.ihe-j.org
--	---	--

IHE Integration Statement

Vendor	Product Name	Version	Date
Karos Health	Rialto for Diagnostic Imaging	2.0	January 19, 2010

This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below:

Profiles Implemented	Actors Implemented	Options Implemented
Registration and Medication History Document Component	Content Consumer	
Sharing Lab Report		View Option
Triage Note		

Internet address for vendor's IHE information: http://www.karoshealth.com/Rialto_Standards.htm		
Links to Standards Conformance Statements for the Implementation		
HL7		
DICOM	http://karoshealth.com/dicom_cs_rialto_20.pdf	
Links to general information on IHE		
In North America: www.ihe.net	In Europe: www.ihe-europe.net	In Japan: www.ihe-j.org