

IHE Integration Statement

Vendor	Product Name	Version	Date
Karos Health	Rialto Vault	1.0	November 01,

This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below:

Profiles Implemented	Actors Implemented	Options Implemented
Access to Radiology Information	Image Manager/Archive	
Audit Trail and Node Authentication	Secure Application	ATNA Encryption
Consistent Presentation of Images	Image Manager/Archive	
Consistent Time	Time Client	
Cross-Enterprise Clinical Documents Share (XDS.b)	Document Repository	
	Document Source	Basic Patient Privacy Enforcement
Cross-Enterprise Document Sharing for Imaging	Document Repository	
	Imaging Document Source	Text Report, PDF Report, Set of DICOM Instances
Cross Enterprise User Authentication	X-Service User	
CT/MR Perfusion with Imaging Contrast	Image Manager/Archive	
Evidence Documents	Image Manager/Archive	
Image Fusion	Image Manager/Archive	
Key Image Note	Image Manager/Archive	
Mammography Image	Image Manager/Archive	
MR Diffusion Imaging	Image Manager/Archive	
NM Image	Image Manager/Archive	
Pathology Workflow	Image Manager/Archive	
Patient Information Reconciliation	Image Manager/Archive	
	Performed Procedure Step Manager	
Radiation Exposure Monitoring	Image Manager/Archive	
Scheduled Workflow	Image Manager/Archive	
	Performed Procedure Step Manager	

Internet address for vendor's IHE information: http://www.karoshealth.com/Rialto_Standards.htm

Links to Standards Conformance Statements for the Implementation

HL7	
DICOM	

Links to general information on IHE

In North America: www.ihe.net	In Europe: www.ihe-europe.net	In Japan: www.ihe-j.org
--	---	--