



## Connect With Your Clients Using Electronic Round-Trips For Orders and Results

### Rialto Consult

Rialto Consult creates an exchange network to improve electronic ordering and results distribution of diagnostic imaging tests between client hospitals, imaging centers and your radiology reading service. Rialto Consult provides a electronic round-trip between a radiology reading service and client institutions.

- Supports secure transmission of images and orders from client institutions to the radiology reading service
- Captures orders and allows client institutions to add notes and other related information, such as scanned documents
- Provides radiologists at the reading service with immediate access to orders, prior imaging results and clinical information, improving turnaround time and reading outcomes
- Keeps physicians updated about the status of their orders
- Distributes results and makes them accessible directly within the client's EMR, RIS or PACS



*Rialto Consult from Karos Health makes it simple for you to extend your reading services to additional hospitals and imaging centres, using your existing RIS/PACS*



## Rialto Consult - A Proven Solution to Extend your Service to more Clients

Rialto Consult is a ready-made solution built on the Rialto clinical information exchange platform, enabling you to grow your radiology reading service through the use of your existing RIS/PACS infrastructure.

### Benefits

#### Efficiency and Cost

- Utilize your existing RIS/PACS infrastructure
- Scale up gracefully with growing data volumes
- Increase provider collaboration

#### Provider Productivity

- Secure access to orders, prior imaging and the complete patient medical history from any connected PACS, RIS or EMR
- Results available are immediately via HL7 integration into native systems
- Track orders and their status

#### Quality of Care

- Diagnose with all relevant information at hand
- Fully comply with privacy legislation and HIPAA
- Fully comply security requirements that does not require VPN

### Rialto Consult

- Implements the IHE Cross-Enterprise Document Sharing (XDS.b) and IHE Cross-Enterprise Document Sharing for Imaging (XDS-I.b) Integration Profiles and includes many actors of the Rialto platform, including a Document Registry, Document Repository, Document Source, Imaging Document Source, Document Consumer, Imaging Document Consumer and Patient Identifier Cross-Reference Manager (PIXV2/PIXV3)
- Ensures that only authorized users can access clinical information and provides complete audit trails in compliance with HIPAA regulation
- Integrates with your existing infrastructure, RIS and PACS
- Provides Single Sign On (SSO) when embedded into existing clinical information systems, augmenting the ability of those systems to display orders, status and clinical information
- Completely web-based, allowing participants access to orders and results from their preferred device
- Promotes collaboration and keeps all participants in the order process up-to-date on order status and results availability through email and other alerts

Learn more about Rialto Consult at [www.karoshealth.com](http://www.karoshealth.com)

**karoshealth**  
*collaboration for life*

7 Father David Bauer Drive  
Suite 201  
Waterloo, ON  
N2L 0A2 Canada

T: +1 519 594 0940  
E: [info@karoshealth.com](mailto:info@karoshealth.com)

[www.karoshealth.com](http://www.karoshealth.com)

### Clinical Information Exchange and Storage Made Simple

Karos Health is focused on elevating the quality of patient care by enabling the sharing and storing of clinical information.

Karos' Rialto platform empowers healthcare enterprises to enable cross-community access to information, facilitating collaboration between healthcare providers and patients.

Rialto is based on open standards that ensure safe and secure handling of patient health information and is backed by a team with decades of experience in healthcare interoperability.